

Weekly Maintenance

The BAQUACIL® CDX™ System includes three key components that work together to save you time and maintenance: BAQUACIL Sanitizer and Algistat, BAQUACIL Oxidizer, and BAQUACIL CDX. See the next page for easy to follow weekly use instructions.

Note: To maintain good water clarity, it is important that the CDX is added IMMEDIATELY following the Oxidizer. Always add BAQUACIL CDX to the skimmer with the pump running.

Important:

- On a weekly basis, it is recommended that you test and adjust your water balance readings to within the ideal ranges discussed on page 8.
- Be sure to use your BAQUACIL Test strips that now test for the presence of Oxidizer.
- Perform the housekeeping tasks discussed on page 14.
- When the water balance levels are within the recommended ranges, follow your BAQUACIL CDX System routine outlined on the next page.

TOP-UP DOSAGE RATES FOR BAQUACIL SANITIZER AND ALGISTAT

	If the level (ppm) is:								
Gallons	40	35	30	25	20	15	10	5	0
5,000	½	¾	¾	1	1¼	1½	1¾	1¾	2
7,500	¾	1	1¼	1½	1¾	2¼	2½	2¾	3
10,000	¾	1¼	1¾	2	2½	2¾	3¼	3¾	4
12,500	1	1½	2	2½	3	3½	4	4½	5
15,000	1¼	1¾	2½	3	3¾	4¼	4¾	5½	6
17,500	1½	2¼	2¾	3½	4¼	5	5¾	6¼	7
20,000	1¾	2½	3¼	4	4¾	5½	6½	7¼	8
22,500	1¾	2¾	3¾	4½	5½	6¼	7¼	8¼	9
25,000	2	3	4	5	6	7	8	9	10
27,500	2¼	3¼	4½	5½	6¾	7¾	8¾	10	11
30,000	2½	3¾	4¾	6	7¼	8½	9¾	10¾	12
32,500	2¾	4	5¼	6½	7¾	9¼	10½	11¾	13
35,000	2¾	4¼	5¾	7	8½	9¾	11¼	12¾	14
37,500	3	4½	6	7½	9	10½	12	13½	15
40,000	3¼	4¾	6½	8	9¾	11¼	12¾	14½	16
42,500	3½	5¼	6¾	8½	10¼	12	13¾	15¼	17
45,000	3¾	5½	7¼	9	10¾	12¾	14½	16¼	18
47,500	3¾	5¾	7¾	9½	11½	13¼	15¼	17¼	19
50,000	4	6	8	10	12	14	16	18	20

All dosages are in pints.